

SEVERAL DIFFERENT CASES

PLACE OF SPECIMEN COLLECTION	PLACE OF SPECIMEN RECEIVING
Ward (bed of patient)	Clinical Pathology Laboratory
Ward (bed of patient)	Microbiological Laboratory
Ward (bed of patient)	Central Clinical Laboratory
Ward (bed of patient)	Forensic Laboratory
Operating Theatre	Clinical Pathology Laboratory
Blood Bank	Haematology Laboratory
Autopsy	Forensic Laboratory
First Aid	Central Laboratory
Satellite Laboratory	Central Laboratory
Day Hospital	Clinical Pathology Laboratory
Private House	Central Laboratory

PLACE OF SPECIMEN COLLECTION	PLACE OF SPECIMEN RECEIVING	SPECIMEN	CONTAINER FOR SPECIMEN	PACKING FOR CONTAINER	REFRIGERATION DURING TRANSIT
Ward (bed of patient)	Pathology Lab (inside the same hospital)	Urine, faeces, sputum, escrementum, blood, sperm	Test tube, cylindrical vessel,	"Biobag", "Bioprotection" "Tran-sbag", "Transeal", "Emo-Express"	Portable refrigerator "Bio-Travel box"
Ward (bed of patient)	Microbiology Lab (inside the same hospital)	Culture, faeces, sputum, escrementum	Test tube, swab, cylindrical vessel	"Biobag", "Bioprotection" "Transeal", "Emo-Express"	Portable refrigerator "Bio-Travel box"
Ward (bed of patient)	Central Clinical Lab (inside the same hospital)	Blood, urine, faeces	Test tube, Vacuum tube, cylindrical vessel	"Biobag", "Bioprotection" "Tran-sbag", "Transeal", "Emo-Express"	Portable refrigerator "Bio-Travel box"
Ward (bed of patient)	Forensic Lab	Blood, urine,	Test tube, Vacuum tube, cylindrical vessel	"Biobag", Bioprotection "Emo-Express" "BioCarrier"	Portable refrigerator "Penguin"
Operating Theatre	Clinical Pathology Lab (inside the same hospital)	Organs, tissues, blood	Test tube, cylindrical vessel	"Biobag", "Bioprotection" "Presto-chiuso"	Portable refrigerator "Penguin" "Bio-Travel box"
Blood Bank	Haematology Lab (inside the same hospital)	Blood	Test tubes, Vacuum test tubes	"Biobag", "Bioprotection" "Emo-Express"	
Autopsy	Forensic Lab (outside hospital)	Organs, tissues, blood	Test tube, cylindrical vessel	"Bioprotection" "Ermetico" "Bio-Carrier"	Portable refrigerator "Penguin"
First Aid	Central Lab (inside the same hospital)	Blood, urine,	Test tube, Vacuum tube, cylindrical vessel	"Biobag", "Bioprotection" "Emo-Express"	
		Organs, tissues, urine,		"Ermetico", "Bio-Carrier"	Portable refrigerator