

Frigostar

LABORATORY REFRIGERATOR

Range of temperature:
from 0 to +12 °C

Automatic defrosting for the
best use of operator's time

High density wall insulation to
reduce the energetic
consumption

Epoxy coated steel structure
with rounded internal corner
for easy cleaning and
disinfection

Inside light activated at
door opening

R404a and CFC-free refrigeration
system



Insulated hinged s/s door
with removable magnetic
rubber gasket to
guarantee a hermetic
closure (triple crystal
door on request)

Forced air to obtain the
best temperature
uniformity and stability
inside the cabinet

Control Panel,
microprocessor
operating for all
functions, acoustic and
visual LED display

Automatic, thermostat
controlled, defrosting

two different control panels RED or SU-X

WHAT TO ORDER

Pannel	Code	Capacity (litres)	Dimensions	Weight (kg)	N° of doors	N° of shelves
RED	SU-X					
89727	89728	370	600x600x1920 h mm	135	1	3
89284	89729	625	720x800x2000 h mm	160	1	3
89377	89730	1365	1440x800x2000 h mm	230	2	6